
## ELITE ELC

### User Instructions

The Sabre Elite range gives low flow oxygen therapy with reduced oxygen consumption.  
12 flow settings give 1 to 6 L/min equivalent oxygen delivery.

#### SAFETY

Follow the instructions of your medical practitioner, the instructions in this leaflet and the instructions for any equipment used with the regulator to ensure optimum patient care and safety.  
DO NOT smoke near the unit at any time as oxygen encourages fire, even in materials which do not normally burn in air.  
DO NOT use oil, wax, polish, or grease on any part of Sabre Elite ELC as they may cause fire or explosion.  
DO NOT use the unit if the 'next service due' date on the waist label has passed.  
Use only MEDICAL OXYGEN with Sabre Elite ELC.  
DO ensure that the cylinder is secure and cannot fall.  
Open and close cylinder valves slowly to avoid oxygen shock.  
ONLY use with a single lumen cannula.  
ALWAYS check that the cannula is unknicked before use.  
ALWAYS close the cylinder valve when you have finished using the Sabre Elite ELC.  
DO NOT tamper with the equipment, as this may endanger patients and will invalidate the warranty.  
DO NOT allow the Sabre Elite ELC to get wet.  
DO NOT cover the breathing hole on the back of the Sabre Elite ELC with labels.  
DO ensure that the gas is vented from the Sabre Elite ELC before removing it from a cylinder.  
If in doubt contact your distributor or Sabre Medical Customer Service Dept. or if you require any guidance in using this product.

#### PROVIDING LOW FLOW OXYGEN THERAPY

Push the nasal cannula over the 'fir-tree' therapy outlet and tug gently to check that it is secure.  
Turn the flow selector knob until the flow setting, stipulated by a medical/nursing practitioner, is displayed in the window by the selector knob.  
Open the cylinder valve slowly and check that the contents gauge shows there is sufficient gas for the therapy to be provided.  
Provide the cannula to the patient and encourage the patient to breathe normally through the nose.  
Monitor the contents gauge regularly during therapy and replace the cylinder when the contents gauge enters the red sector.  
Ensure that the cannula remains in place throughout therapy. If it is displaced no therapy will be provided. When a patient inhales, Sabre Elite ELC can be heard providing pulses of oxygen. During exhalation, or if the cannula is displaced from the nose the unit is silent and there is no oxygen flow.  
Note: The best oxygen savings will be obtained if there is one pulse per breath. It is worth experimenting with positioning of the cannula in the nose to achieve this.

#### AFTER USE

Remove the nasal cannula from the patient, check the amount of gas remaining in the cylinder and close the cylinder valve. Re-fill empty cylinders.

#### CHARGING AN ELITE ELC CYLINDER

Follow the instructions provided with the Sabre Charging Adaptor to connect to the charging cylinder and charge the Sabre Elite ELC. Replace any dust caps and store the charging adaptor in a sealed plastic bag or container.

#### CLEANING

DO NOT dismantle the unit.  
DO NOT allow moisture to enter the unit.  
DO NOT use polish including spray polish, oil, or any other cleaning agent or lubricant.  
Dust the exterior surfaces with a clean dry lint-free cloth.  
A clean, slightly moist, lint free cloth may be used to remove persistent marks.

DO NOT pack the equipment until completely dry.

#### SERVICE REQUIREMENTS

The unit must be serviced by the date marked on the waist label.  
Sabre Elite ELC must be serviced by Sabre Medical.  
Cylinders on direct units must be re-tested every 5 years in accordance with gas transport regulations. The unit will be serviced at the same time.

#### PRINCIPLE OF OPERATION

Oxygen from the cylinder is reduced to a constant pressure and used to fill a reservoir at a rate set by the flow setting knob.  
During exhalation, there is no oxygen flow out of the unit, so the reservoir fills up to a level determined by the flow setting.  
As an inhalation begins, the small suction at the nose makes a valve open. A timed pulse of oxygen is delivered from the reservoir to the patient. The timing is designed so inhaled gas flows deep into the lungs.  
Flows are set so that the oxygen delivered into the lung will be approximately equivalent to the same constant flow.  
Flows into the reservoir are about a third of the flow setting. This, combined with the number of settings will normally give increases in duration of 3 to 5.

#### TECHNICAL DESCRIPTION

Gas supply: 10 to 200 bar  
Therapy outlet: Universal 'fir-tree' type.  
Therapy flow accuracy: +/- 10%  
Outlet pressure: 4.2 bar (61 PSI) nominal.  
Operating and storage limits: -20°C to +60°C, humidity up to 95% RH non condensing.

#### SPARES AND ACCESSORIES

Nasal cannula: 1025407  
Charging adaptors for refilling PORTAFLOW from larger cylinders:  
UK fitting: 2003558  
French fitting (NFE-29-650): 1018432  
US fitting (CGA 540): 1018433  
German fitting (DIN477 type 9): 1018434

#### SABRE MEDICAL PRODUCTS

Sabre Medical Products offer a wide range of oxygen therapy, resuscitation and medical gas equipment. For more details contact Sabre Customer Service

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