

# Liquid Oxygen Use Time Chart: Nutzungsdauer

Model	0H	0.12	0.25	0.50	0.75	1.00	1.50	2.00	2.50	3.00	3.50	4.00	5.00	6.00	8.00	10.00	12.00	15.00
Model	0H	0.12	0.25	0.50	0.75	1.00	1.50	2.00	2.50	3.00	3.50	4.00	5.00	6.00	8.00	10.00	12.00	15.00
HELIOS Plus	21.9	21.9	17.8	12.0	7.6	18.4	12.2	9.2	7.3	6.1	5.2	4.6	-	-	-	-	-	-
HELIOS Marathon	48.2	-	-	-	-	12.1	29.3	6.1 22.0	17.6	4.0 14.6	-	3.0 11.0	2.4	2.0	-	-	-	-
Spirit 300	15.8	-	-	-	-	15.1	10.7	2.4 8.0	-	5.4	-	4.0	3.2	-	-	-	-	-
Spirit 600	31.7	-	-	-	-	30.3	21.4	4.8 16.1	-	10.7	-	8.0	6.4	-	-	-	-	-
Spirit 1200	60.0	-	-	-	-	57.4	40.5	9.2 30.4	-	20.3	-	15.2	12.2	-	-	-	-	-
Sprint	31.7	-	18.0	12.2	9.3	7.5	5.5	4.3	3.6	3.1	-	2.4	2.0	1.7	-	-	-	-
Stroller	60.0	-	34.0	23.0	17.4	14.1	10.2	8.0	6.6	5.6	-	4.4	3.6	3.0	-	-	-	-
Hi Flow Stroller	60.0	-	-	23.0	-	14.1	-	8.0	6.6	5.6	-	4.4	-	3.0	2.3	1.9	1.6	1.4
C1000	54.5	-	30.8	21.5	16.5	13.4	9.7	7.7	6.3	5.4	-	4.1	3.3	2.8	-	-	-	-
C1000T	54.5	-	-	21.5	-	13.4	9.7	7.7	-	5.4	-	4.1	3.3	2.8	2.1	1.7	-	1.2
C500	27.9	-	27.9	18.3	12.2	9.1	6.1	4.6	3.7	3.0	-	2.3	1.8	1.5	-	-	-	-
C550	27.9	-	27.9	18.3	12.2	9.1 12.1	6.1 8.1	4.6 6.1	3.7 4.9	3.0 4.0	-	2.3 3.0	1.8 2.4	1.5 2.0	-	-	-	-
Reservoirs																		
Liberator 10	21.4	-	19.8	13.3	8.7	6.5	4.3	3.2	2.6	2.2	-	1.6	1.3	1.1	-	-	-	-
Liberator 20	34.6	-	34.6	26.9	17.5	13.0	8.7	6.5	5.2	4.3	-	3.3	2.6	2.2	1.6	1.3	1.1	0.9
Liberator 30	49.9	-	49.9	38.9	25.3	18.8	12.6	9.4	7.5	6.3	-	4.7	3.8	3.1	2.4	1.9	1.6	1.3
Liberator 37	61.3	-	61.3	47.7	31.1	23.1	15.4	11.6	9.3	7.7	-	5.8	4.6	3.9	2.9	2.3	1.9	1.5
Liberator 45	74.6	-	74.6	58.0	37.8	28.1	18.8	14.1	11.3	9.4	-	7.0	5.6	4.7	3.5	2.8	2.3	1.9
Liberator 60	90.0	-	90.0	74.7	48.7	36.2	24.2	18.1	14.5	12.1	-	9.1	7.2	6.0	4.5	3.6	3.0	2.4
Companion 21	26.5	-	26.5	22.6	17.6	12.9	8.8	6.5	5.1	4.3	-	3.2	2.5	2.1	-	-	-	-
Companion 31	39.2	-	39.2	33.4	26.1	19.1	12.9	9.6	7.6	6.3	-	4.7	3.7	3.1	-	-	-	-
Companion 41	51.8	-	51.8	44.2	34.5	25.2	17.1	12.7	10.0	8.3	-	6.2	5.0	4.1	-	-	-	-
HELIOS 36	72.0	-	-	44.8	29.2	21.7	14.5	10.9	-	7.2	-	5.4	4.3	3.6	2.7	2.2	-	-
HELIOS 46	92.0	-	-	57.3	37.3	27.8	18.5	13.9	-	9.3	-	6.9	5.6	4.6	3.5	2.8	-	-
Universal 36	72.0	-	-	44.8	29.2	21.7	14.5	10.9	-	7.2	-	5.4	4.3	3.6	2.7	2.2	-	-
Universal 46	92.0	-	-	57.3	37.3	27.8	18.5	13.9	-	9.3	-	6.9	5.6	4.6	3.5	2.8	-	-

Hinweis:  
Die Nutzungsdauer der tragbaren Einheiten sind in Stunden dargestellt.  
Die Nutzungsdauer der stationären Einheiten sind in Tagen dargestellt.  
Die Nutzungsdauer der Dosiersysteme basieren auf 20 bpm (Atemzüge pro Minute) und sind.